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MHC multimer assays: Optimized conditions allowing enumeration of low-affinity, antigen specific CD4+ and CD8+ T cells by flow cytometry

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
Reliable monitoring of cellular immune responses is important in vaccine and immunotherapeutic development. Flow analysis using conventional MHC multimers like Tetramers and Pentamers has made a great impact on this field, enabling enumeration and phenotypic characterization of antigen-specific T cells. However, it is often difficult to obtain well-separated, distinct negative and positive cell populations when attempting to detect T cells of low affinity for the MHC multimer.

The use of MHC Dextramers can overcome these problems, and provide reliable detection of low-affinity CD8+, CD4+ and iNKT cells. MHC Dextramers are an improved form of MHC multimers that improves staining intensity, resolution and signal-to-noise dramatically under these circumstances.

Cancer vaccines typically stimulate T cells with low affinity for the cognate pMHC class I complex. We show that MHC Dextramers, in contrast to Tetramers, can efficiently stain low-affinity T cells. Thus, using Dextramers, staining is achieved for affinities as low as 500 μ M, and effective separation of the negative and positive cell populations is achieved for affinities as low as 250 μ M.

pMHC class II complexes typically have low affinity for their cognate TCRs, making it difficult to enumerate CD4+ cells involved in e.g. autoimmune or inflammatory disease using conventional MHC multimers. We show that class II Dextramers allow detection of CD4+ T cells that cannot be detected using tetramers.

Finally, CD1d complexes (MHC-like complexes comprising glycolipids rather than peptides) often have low affinity for the TCR of iNKT cells. The use of CD1 Dextramers may be used to enumerate these cells.

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Immudex
IMMUNOLOGY TAKEN TO A HIGHER LEVEL

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