

## Instructions for Presenters in “One Minute Presentation”

The abstracts which were only adopted to poster presentation, and categorized into program number “P”, and are entitled to do the “One Minute Presentation” orally in the relevant workshop.

Ex.1-Q-W99-99-**P** (The mark “O/P” means Oral and Poster presentation and is not subjected to the “One Minute Presentation”.)

**50 seconds** are allotted for each presentation and presenters may present his/her work briefly. The PC slide for the presentation is limited to 1 slide.

### 【How to prepare the PC slide for “One Minute Presentation”】

1. Prepare one horizontal A4 slide, in PDF format with Adobe Acrobat V.7 or 8 or 9.
2. If you have any trouble and difficulties in formatting the file to PDF, please contact to the Support Desk in Hibino Media Technical ( [jsi40-omp@media-t.co.jp](mailto:jsi40-omp@media-t.co.jp)).
3. No movies, no animation and no sound effects are permitted.
4. The slide should be followed as given below format, in order to look clearly from audience.
5. Font: Mac users or PC users other than Windows must embed any fonts that are not supported by Windows in the PDF.

1-A-W1-1-P Your Name (Meneki Univ Med Sch)	
<b>XXXXX is essential for YYYYY in ZZZZZ</b>	← Title of your presentation
XXXXX is known for its involvement in AAAAA inBBBBB (Nature 1950). However, the role of XXXXX in ZZZZZ is unknown.	
← Brief Background	
Here we show that	← Key results
1) The expression of XXXXX in PPP cells is up-regulated in ZZZZZ patients; 2) Genetic ablation of XXXXX in mice reduces the YYYYY response of PPP cells; and 3) Neutralizing antibodies to XXXXX block YYYYY response of human PPP cells.	
<div style="border: 2px solid green; padding: 10px; text-align: center;"><h2><u>ONE OR TWO KEY DATA</u></h2></div>	
Conclusions: ← Brief conclusions (optional)	
- XXXXX is essential for YYYYY response of PPP cells in mouse and human. - XXXXX may be a good target for diagnosis and therapy of ZZZZZ.	

6. Please note that the name of the PDF file will be registered automatically to the Program Number.